Varioskan LUX軟體操作說明 v4.1

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點選Report方塊(圖1-7)
可選擇列印與報告輸出格式(Table, List)
Open in excel (excel檔)
Export report to file (Text 檔)
Print

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Action	Description
Absorbance	Measures absorbance.
Absorbance Spectrum	Measures the absorbance spectrum for a wavelength range.
Fluorescence	Measures fluorescence intensity.
Fluorescence Spectrum	Measures the fluorescence spectrum for a wavelength range.
Luminescence	Measures luminescence.
Luminescence Spectrum	Measures the luminescence spectrum for a wavelength range.
AlphaScreen*	Measures the AlphaScreen* signal.
TRF	Measures time-resolved fluorescence.
TRF Spectrum	Measures time-resolved fluorescence spectrum for a wavelength range.
Kinetic Loop	Executes sub-steps several times in defined time intervals in a kinetic measurement.
Well Loop	Executes sub-steps for as many wells at a time as you have selected as a well count.
Area Selection	Executes sub-steps for only part of the defined wells in the plate layout. Area definition is not necessary when all of the defined wells in the plate layout are measured.
Shake	Shakes the microplate to mix the liquid in the wells.
Dispense	Dispenses a given volume of liquid into wells.
Pause	Pauses the protocol.
Run Plate Out / In	Runs the plate in or out in the middle of a protocol.

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
Example: To add a series of standard samples (concentrations 5, 10, 50, 100 and 500 ug/mL), with two replicates side by side, select the pipette content parameters as shown below, and add samples to the plate by painting the wells with the virtual pipette:

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