CFX Real Time PCR system

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
The 'Plate Loading Guide' dialog box is shown with the following fields and annotations:

- Select Fluorophores...**: 選取實驗使用的螢光試劑
- Sample Type**: Unknown (選取樣品種類(unknown, standard, NTC, PC, NC, NRT))
- Target Names**:
 - Load ☒ FAM <none> +
 - Load ☒ HEX <none> +
 - Load ☒ Texas <none> +
 - Load ☒ Cy5 <none> +
 - Load ☒ Quasar <none> +
 選取此次實驗所用的螢光試劑，並註記 target name 名稱
- Sample Names**: Load ☐ <none> + 輸入樣品名稱
- Biological Group**: Load ☐ <none> + 輸入組別名稱
- Replicate #**: Load ☐ 1 輸入重複性樣品設定(Horizontal 水平重複/Vertical 垂直重複)
- Technical Replicates**: Show Technical Replicates (info icon)
- Experiment Settings...**: 設定 ref gene 與 control sample (圖十)
- Clear Replicate #**: 清除所選的重複
- Clear Wells**: 清除所選的輸入

Buttons: OK, Cancel

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選擇 Reference gene

選擇 Control Sample

The image shows two side-by-side 'Experiment Settings' windows. The left window is titled '選擇 Reference gene' and the right window is titled '選擇 Control Sample'. Both windows have a 'Targets' tab and a 'Samples' tab. The 'Samples' tab is active in both. The left window has a table with columns: Name, Full Name, Reference, and Select To Remove. The right window has a table with columns: Name, Full Name, Control, and Select To Remove. In the left window, 'GAPDH' is selected in the 'Reference' column. In the right window, sample '1' is selected in the 'Control' column. Both windows have a 'New:' field, an 'Add' button, and a 'Remove checked item(s)' button. There are also checkboxes for 'Show Analysis Settings' and 'Exclude the following sample types from Gene Expression analysis:'. The 'Exclude' section has checkboxes for NTC, NRT, Negative Control, Positive Control, and Standard.

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The image shows a 'Starting Concentration' dialog box. It has the following fields and controls:

- Starting Concentration: 1.00E+06 (with a red arrow pointing to it from the label '標準品的起始濃度')
- Replicates from: 1 (with a red arrow pointing to it from the label '序列稀釋次數')
- to: 3 (with a red arrow pointing to it from the label '序列稀釋次數')
- Dilution Factor: 10.000 (with a red arrow pointing to it from the label '稀釋倍數')
- Radio buttons: Increasing (unselected) and Decreasing (selected)
- A dropdown menu: <All>
- Buttons: Cancel and Apply

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