BIO-RAD
CFX Opus

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
The MIQE Guidelines:
Minimum Information for Publication of Quantitative
Real-Time PCR Experiments

Stephen A. Bustin,^{1*} Vladimir Benes,² Jeremy A. Garson,^{3,4} Jan Hellemans,⁵ Jim Huggett,⁶
 Mikael Kubista,^{7,8} Reinhold Mueller,⁹ Tania Nolan,¹⁰ Michael W. Pfaffl,¹¹ Gregory L. Shipley,¹²
 Jo Vandesompele,⁵ and Carl T. Wittwer^{13,14}

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MIQE 準則

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
MIQE 準則

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MIQE 準則

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1

檢核項目	重要性	說明
EXPERIMENTAL DESIGN		
Definition of experimental and control groups	E	敘述實驗設計組別、樣本數以及執行的實驗室等
Number within each group	E	
Assay carried out by core lab or investigator's lab?	D	
Acknowledgement of authors' contributions	D	
SAMPLE		
Description	E	敘述樣本來源(例如組織或細胞)、細胞數(量)、處理方式(多少時間內完成)、保存環境(溫度)
Volume/mass of sample processed	D	
Microdissection or macrodissection	E	
Processing procedure	E	
If frozen - how and how quickly?	E	
If fixed - with what, how quickly?	E	
Sample storage conditions and duration (especially for FFPE samples)	E	

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MIQE 準則

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qPCR TARGET INFORMATION		
If multiplex, efficiency and LOD of each assay.	E	於同一試管做多基因定量才需要確定該項目
Sequence accession number	E	先以資料庫或軟體確定基因序列、比對、二級結構分析、Primer 對應於 exon 位置等資訊
Location of amplicon	D	
Amplicon length	E	
<i>In silico</i> specificity screen (BLAST, etc)	E	
Pseudogenes, retropseudogenes or other homologs?	D	
Sequence alignment	D	
Secondary structure analysis of amplicon	D	
Location of each primer by exon or intron (if applicable)	E	
What splice variants are targeted?	E	

qPCR OLIGONUCLEOTIDES		
Primer sequences	E	Primer 或 Probe 合成的序列、廠商、純化方法等
RTPrimerDB Identification Number	D	
Probe sequences	D ³	
Location and identity of any modifications	E	
Manufacturer of oligonucleotides	D	
Purification method	D	

qPCR PROTOCOL		
Complete reaction conditions	E	敘述定量 PCR 實驗反應條件
Reaction volume and amount of cDNA/DNA	E	cDNA 量與反應體積
Primer, (probe), Mg ⁺⁺ and dNTP concentrations	E	Primer、鎂離子濃度、聚合酶濃度、緩衝液組成、其他添加物等
Polymerase identity and concentration	E	
Buffer/kit identity and manufacturer	E	
Exact chemical constitution of the buffer	D	
Additives (SYBR Green I, DMSO, etc.)	E	
Manufacturer of plates/tubes and catalog number	D	分析盤(管)廠牌和型號
Complete thermocycling parameters	E	反應條件(溫度/時間/cycles)
Reaction setup (manual/robotic)	D	反應設定
Manufacturer of qPCR instrument	E	儀器廠牌機型

qPCR VALIDATION		
Evidence of optimisation (from gradients)	D	是否有測試 annealing 溫度梯度
Specificity (gel, sequence, melt, or digest)	E	確認反應專一性
For SYBR Green I, Cq of the NTC	E	是否 NTC 控制組有產物
Standard curves with slope and y-intercept	E	建立標準曲線 ⁴ :由斜率推算 PCR 效率、計算信賴區間與標準差以及 r ²
PCR efficiency calculated from slope	E	
Confidence interval for PCR efficiency or standard error	D	
r ² of standard curve	E	
Linear dynamic range	E	求得線性關係範圍 ⁵ : 計算信賴區間、最低極限時的 Cq 變異數以及是否確認 LOD 偵測極限
Cq variation at lower limit	E	
Confidence intervals throughout range	D	
Evidence for limit of detection	E	
If multiplex, efficiency and LOD of each assay.	E	於同一試管做多基因定量計算效率和偵測極限

DATA ANALYSIS		
qPCR analysis program (source, version)	E	分析程式/工具(版本)、分析方法 (threshold 或 regression)以及敘述樣本盤排列方式等
Cq method determination	E	
Outlier identification and disposition	E	
Results of NTCs	E	NTC 控制組是否有產物
Justification of number and choice of reference genes	E	Reference gene 的數目和選擇
Description of normalisation method	E	敘述標準化方法
Number and concordance of biological replicates	D	生物樣本數
Number and stage (RT or qPCR) of technical replicates	E	單一生物樣本在儀器偵測的 n 值
Repeatability (intra-assay variation)	E	計算重複性、再現性、使用何種統計方法、統計軟體(版本)
Reproducibility (inter-assay variation, %CV)	D	
Power analysis	D	
Statistical methods for result significance	E	
Software (source, version)	E	
Cq or raw data submission using RDML	D	以 RDML 格式提供原始數據

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
MIQE 準則

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Biological Replicates

Technical Replicates


Fig. 6. Experimental replicates. All experiments should be designed with a combination of biological and technical replicates. This illustrates a simple experiment with triplicate biological samples from control (■) and treatment/experimental (■) conditions. For each biological sample, three technical replicates are recommended for the gene of interest as well as for the reference gene(s). This results in a total of at least 36 samples plus the duplicate NTC (■) for a total of >40 wells.

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MIQE 準則

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Sample	Assay
Bile salts	Cholic acid Deoxycholic acid
Polysaccharides	Polygalacturonic acid
Lipids	Cholesterol hydrochloride
Algae	Alginic acid
Clay	Montmorillonite
Blood	Heparin Hematin EDTA Serum
Soil	Humic acid
Textile	Indigo
Wine	Tannic acid
Plants	Cellulose Pectin (for fiber control)
Hair, tissues	Melanin
Bones, teeth	CaCl ₂

Sample	Assay
Cell lysates	MEM+FBS PBS Trypsin
Sample preparation	NaAc NaCl EtOH Isopropanol TRIzol InstaGene Matrix
Miscellaneous	Green tea Chocolate Dust SDS DMSO DTT
Body fluids	Spermidine Urea

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Genmall Biotechnology Co., Ltd.