# Constant Systems Ltd

since 1989

## Continuous High Pressure Cell Disrupter

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
Algae, animal cells, bacteria, environmental samples, fungi, yeast, mammalian cells and tissue, parasites, plant cells and tissues and viruses.


- Protein & DNA Extraction
  - DNA Manipulation
  - Enzyme Release
  - Single Cell Isolation
  - Selective Breakage
  - Tissue Disruption
- 能獲得大量產物
  - 避免在破碎過程中變性
  - 避免有害物質產生而影響後續實驗

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

Candida Albicans Yeast  
Form Unbroken


Candida Albicans Yeast  
Form Broken

# TS 0.75KW 系統

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*Graph showing a typical controlled pressure cycle in the cell disrupter*

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

# Common Cell Disruption Models


## Yeast - *Saccharomyces Cerevisiae*

**Protein Release** - more than 99% after one pass at 40 Kpsi

**Membrane Breakage at ATPase activity** - best at 30 Kpsi

**Protein Purification** - best at 20 Kpsi

**Recombinant Protein** - best at 13.5 Kpsi

# Common Cell Disruption Models

## Bacteria - *Escherichia Coli* (E-Coli)


<b>Membrane Protein</b>	- best between 27 and 35 Kpsi
<b>Cytoplasmic Protein</b>	- best between 20 and 27 Kpsi
<b>Inclusion Bodies</b>	- best between 15 and 20 Kpsi
<b>Release of DNA</b>	- best below 15 Kpsi

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
	Flow	Pressure	Flow Rate	Features
<b>One Shot Model 0.75KW</b>	One Shot	Up to 40,000psi	0.5ml - 20ml per shot	Research
<b>TS 0.75kW</b>	Continuous Flowing	Up to 40,000psi	40ml - 115ml /min	Research
<b>TS 1.1KW</b>	Continuous Flowing	Up to 40,000psi	100ml - 255ml /min	Pilot Scale
<b>TS 2.2KW</b>	Continuous Flowing	Up to 40,000psi	190ml - 310ml /min	Pilot/Production Scale
<b>TS 4KW</b>	Continuous Flowing	Up to 40,000psi	405ml - 565ml /min	Production Scale

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

# Examples of High Pressure Disruption

## 2. The effect of pressure on the cell disruption of *Saccharomyces Cerevisiae* expressing protein X

*Completed at an independent and confidential customer site, 2003.*

### Preparation:

Cell pellet was resuspended in chilled 20mM Tris/1 mM EDTA, pH 7.98. The required lysis buffer was calculated as follows: (46.82% w/v/20%w/v) x 233ml = 545ml. The resuspended cells were mixed for 30 minutes at 6°C. pH and conductivity of resuspended paste; 6.98 at 12.2°C and 1.448mS/cm at 13.8°C. OD<sub>600</sub> of resuspended cells; 76.

### Method 1:

40ML of the resuspended cells were passed through a Constant Systems disrupter once at the following pressures; 8, 10, 15, 20, 25 and 27kpsi.

Pressure	Protein mg/ml
Non Disrupted	0.88
8	2.20
10	2.89
15	4.16
20	4.34
25	4.41
27	4.61

## Result after 1 pass

Sample	% Cell Disruption
Non Disrupted	0
8	52%
10	68%
15	73%
20	87%
25	96%
27	98.9%

# Examples of High Pressure Disruption


## 1. The effect of pressure on the cell disruption of *Esherichia coli* expressing soluble protein X

*Completed at an independent and confidential customer site, 2003.*

### Preparation:

Cell paste (87g WCW) resuspended in 870ml 20mM TrisHCl/5 mM EDTA, pH 8.0. 40ml of the suspension was passed through Constant Cell Disruption Systems (Z Plus Series) once with the following disruption pressures: 10, 17, 20, 24, 27, 30, 33 and 35kPSI respectively. Samples were centrifuged at 8000 rpm for 20 minutes at 4°C. Determined protein released with coomassie assay. Thereafter repeat passes of the same sample were passed through the cell disrupter to determine if multiple passes improved efficiency in protein yield between passes.

### Results 1 – Protein Release:

*For this particular strain of e.coli 27kpsi was the optimal pressure.*

Pressure (kpsi)	Protein (ug/ml)
0	414
10	3221
17	3838
20	4106
24	4162
27	4261
30	4247
33	4240
35	4236

Number of Passes	Protein ug/ml
1	4261
2	4840
3	4884
4	4488

Cell Disruption following 1 pass through the cell disrupter = 99.87%.  
99.87% breakage would leave 0.13% intact =  $7.2 \times 10^7$  cfu/ml

Cell Disruption following 2 pass through cell disrupter = 99.997%.  
99.997% breakage would leave 0.003% intact =  $1.7 \times 10^6$  cfu/ml.

# Constant Systems 使用客戶

- Pfizer -UK
- NIH -USA
- University of Cambridge -UK
- Tokyo University -Japan
- SmithKline Beecham -Belgium
- Academia Sinica -Taiwan
- NHRI-Taiwan
- And More...

Thanks for Your Attending.