

下載 QuantStudio™ 3 即時聚合酶鏈鎖反應儀軟體方式

1. 連結至[下載網頁](#)，點選 Get software
QuantStudio 3/5 Real-Time PCR Software



The QuantStudio 3 and 5 real-time PCR system workflow is supported by the following software applications.

Design & Analysis 1 (DA1) software and Touchscreen Instrument Operating software

This desktop instrument operating software enables users to operate QuantStudio 3 and 5 instruments. Customers who need to operate the QuantStudio 3/5 system within the desktop software or to access the Securities, Audits, and e-Signature module for the QuantStudio 5 system must use the Design and Analysis v1.x series.

[Get software](#)

2. 因原廠規定須登入帳號後下載，如您尚無帳號請先完成註冊後再登入下載。

ThermoFisher
SCIENTIFIC

登入您的帳戶

電子郵件/使用者名稱

使用者名稱為必填欄位

下一步

無主登入嗎？

沒有帳戶嗎？

建立帳戶

3. 登入後於左方頁籤找到軟體下載點，點選右上 Download 完成檔案下載，並依照提示完成安裝即可。

DataConnect

New Experience (Beta)

Manage and share your files and projects

Create a group Upload files New folder Create project

Previous versions Download Move Delete Share

File Name	File Type	Run Date	Modified D...
DA 1.5.3 Installer (Win 10).zip	Shared ZIP		04/Nov/2023
DA 1.5.3 Release Notes.pdf	Shared PDF		11/Oct/2023

DA 1.5.3 Installer (Win 10).zip

Data

Size: 409.6 MB

Updated by: techsupport@thermofisher.com

Last modified: 04/Nov/2023 02:10:02 AM

Run date:

Data type:

File type: ZIP

Owner: QuantStudio 1/3/5

Custom Information

eulaUrl: https://apps.thermofisher.com/dhap-resources/eula/DAD-qPCR-terms.html

New label Add description