

# 下載 QuantStudio™ 3/5 即時聚合酶鏈鎖反應儀軟體方式

1. 連結至[下載網頁](#)，點選 Get software

QuantStudio 3/5 Real-Time PCR Software



The QuantStudio 3 and 5 real-time PCR system workflow is supported by the following software applications.

#### Design & Analysis 1 (DA1) software and Touchscreen Instrument Operating software

This desktop instrument operating software enables users to operate QuantStudio 3 and 5 instruments. Customers who need to operate the QuantStudio 3/5 system within the desktop software or to access the Securities, Audits, and e-Signature module for the QuantStudio 5 system must use the Design and Analysis v1.x series.

[Get software](#)

下載網頁連結



2. 因原廠規定須登入帳號後下載，如您尚無帳號請先完成註冊後再登入下載。

3. 登入後於左方頁籤找到軟體下載點，在 DA 1.6.1 Installer(Win 10)按右鍵，出現 Download 下載檔案，於電腦上安裝即可。